



Screw guide



EXTERNAL HEX UNIVERSAL

compatible with
Nobel Brånemark®

SOCKET
UniGrip®



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.027	Screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.028	Screw UG RP	2,53	0°	7,35	1,42	3,12	2,81	M2,0	35 Ncm
	19.029	Screw UG RP	2,85	0°	7,10	1,40	2,70	3,00	M2,5	35 Ncm
	19.227	DLC screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.228	DLC screw UG RP	2,53	0°	7,35	1,42	3,12	2,81	M2,0	35 Ncm
	19.229	DLC screw UG RP	2,85	0°	7,10	1,40	2,70	3,00	M2,5	35 Ncm
	19.627	TIN screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.628	TIN screw UG RP	2,53	0°	7,35	1,42	3,12	2,81	M2,0	35 Ncm
	19.629	TIN screw UG RP	2,85	0°	7,10	1,40	2,70	3,00	M2,5	35 Ncm




TRI-LOBE

compatible with
NobelReplace Select™

SOCKET
UniGrip®



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.004	Screw UG NP	2,45	0°	8,30	1,70	2,60	2,75	M1,6	35 Ncm
	19.005	Screw UG RP & WP	2,45	0°	9,90	1,50	3,12	2,81	M2,0	35 Ncm
	19.204	DLC screw UG NP	2,45	0°	8,30	1,70	2,70	3,00	M2,5	35 Ncm
	19.205	DLC screw UG RP & WP	2,45	0°	9,90	1,50	2,60	2,75	M1,6	35 Ncm
	19.604	TIN screw UG NP	2,45	0°	8,30	1,70	3,12	2,81	M2,0	35 Ncm
	19.605	TIN screw UG RP & WP	2,45	0°	9,90	1,50	2,70	3,00	M2,5	35 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



ACTIVE HEX

compatible with
NobelActive® / Replace® CC

SOCKET
UniGrip®



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.027	Screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.030	Screw UG RP	2,53	54°	7,35	1,57	2,97	2,81	M2,0	35 Ncm
	19.068	Screw UG 3.0	1,82	55°	6,72	1,82	2,45	2,45	M2,5	15 Ncm
	19.227	DLC screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.230	DLC screw UG RP	2,53	54°	7,35	1,57	2,97	2,81	M2,0	35 Ncm
	19.268	DLC screw UG 3.0	1,82	55°	6,72	1,82	2,45	2,45	M2,5	15 Ncm
	19.627	TIN screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 Ncm
	19.630	TIN screw UG RP	2,53	54°	7,35	1,57	2,97	2,81	M2,0	35 Ncm
	19.668	TIN screw UG 3.0	1,82	55°	6,72	1,82	2,45	2,45	M2,5	15 Ncm




OCTAGON

compatible with
Straumann® Tissue level & synOcta®

SOCKET
Torx 06



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.009	Screw Torx on abutment	2,50	45°	4,50	1,20	0	3,30	M2,0	15 Ncm
	19.026	Screw Torx	2,58	75°	6,70	3,20	1,68	2,00	M2,0	35 Ncm
	19.070	Screw Torx NNC	2,10	75°	7,10	1,85	3,35	1,70	M1,6	35 Ncm
	19.609	TIN screw Torx on abutment	2,50	45°	4,50	1,20	0	3,30	M2,0	15 Ncm
	19.626	TIN screw Torx	2,58	75°	6,70	3,20	1,68	2,00	M2,0	35 Ncm
	19.670	TIN screw Torx NNC	2,10	75°	7,10	1,85	3,35	1,70	M1,6	35 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



CONICAL BL

compatible with
Straumann® Bone level

SOCKET
Torx 06



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.032	Screw Torx NP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.033	Screw Torx RP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.037	Screw Torx NC/RC	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.098	Screw for SRA	2,20	55°	3,65	1,65	0,20	1,70	M1,4	15 Ncm
	19.632	TIN screw Torx NP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.633	TIN screw Torx RP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.637	TIN screw Torx NC/RC	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 Ncm
	19.698	TIN screw for SRA	2,20	55°	3,65	1,65	0,20	1,70	M1,4	15 Ncm



CONICAL BLX

compatible with
Straumann® BLX

SOCKET
Torx 06



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.054	Screw Torx	2,10	90°	6,10	1,02	4,26	1,84	M1,6	35 Ncm
	19.700	TIN screw Torx	2,10	90°	6,10	1,02	4,26	1,84	M1,6	35 Ncm




NEO GM

compatible with
Neodent® Grand Morse

SOCKET
GM



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.575	Screw NEO	1,86	60°	9,75	1,50	7,65	2,15	M1,6	20 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



INTERNAL HEX CONIC

compatible with
Astra Tech Osseospeed™

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.015	Screw hex NP (3.0)	2,10	30°	7,50	2,00	3,00	2,50	M1,4	15 Ncm
	19.016	Screw hex RP	2,35	0°	8,25	2,00	3,22	3,02	M1,6	20 Ncm
	19.017	Screw hex WP	2,35	0°	8,30	2,00	3,80	3,50	M2,0	25 Ncm
	19.018	Screw on abutment	2,30	0°	4,80	1,90	0,00	2,90	M1,4	15 Ncm
	19.615	TIN screw hex NP (3.0)	2,10	30°	7,50	2,00	3,00	2,50	M1,4	15 Ncm
	19.616	TIN screw hex RP	2,35	0°	8,25	2,00	3,22	3,02	M1,6	20 Ncm
	19.617	TIN screw hex WP	2,35	0°	8,30	2,00	4,80	3,50	M2,0	25 Ncm
	19.618	TIN screw on abutment hex	2,30	0°	4,80	1,90	0,00	2,90	M1,4	15 Ncm




CONIC EVO

compatible with
Astra Tech implant system™ EV

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.058	Screw hex on Uniabutment EV	2,35	18°	4,49	1,71	0,00	2,50	M1,8	15 Ncm
	19.059	Screw hex EV 3.6	2,15	60°	7,46	2,12	1,34	4,00	M1,6	25 Ncm
	19.060	Screw hex EV 4.2	2,35	40°	7,64	1,85	1,79	4,00	M1,8	25 Ncm
	19.061	Screw hex EV 4.8/5.4	2,35	0°	7,60	1,60	2,00	4,00	M2,0	25 Ncm
	19.071	Screw hex EV 3.0	2,15	52°	7,68	1,60	3,83	3,85	M1,4	25 Ncm
	19.658	TIN screw hex on Uniabutment EV	2,35	18°	4,49	1,71	0,00	2,50	M1,8	15 Ncm
	19.659	TIN screw hex EV 3.6	2,15	60°	7,46	2,12	1,34	4,00	M1,6	25 Ncm
	19.660	TIN screw hex EV 4.2	2,35	40°	7,64	1,85	1,79	4,00	M1,8	25 Ncm
	19.661	TIN Screw hex EV 4.8/5.4	2,35	0°	7,60	1,60	2,00	4,00	M2,0	25 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



INTERNAL HEX FD

compatible with
Friadent Xive®

SOCKET
HEX 1,22mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.024	Screw hex	2,10	45°	9,80	3,06	3,29	3,45	M1,6	25 Ncm



19.624	TIN screw hex	2,10	45°	9,80	3,06	3,29	3,45	M1,6	25 Ncm
--------	---------------	------	-----	------	------	------	------	------	--------



INTERNAL ANK

compatible with
Dentsply Ankylos® C/X

SOCKET
HEX 1,00mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.048	Screw hex	2,35	60°	7,87	2,37	3,40	3,10	M1,8	15 Ncm



19.050	Screw hex	2,35	60°	9,37	2,02	3,35	3,40	M1,8	15 Ncm
--------	-----------	------	-----	------	------	------	------	------	--------



19.648	TIN screw hex	2,35	60°	7,87	2,37	2,40	3,10	M1,8	15 Ncm
--------	---------------	------	-----	------	------	------	------	------	--------



EXTERNAL HEX USA

compatible with
3i® Osseotite®

SOCKET
HEX 1,20mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.034	Screw hex	2,55	0°	7,85	1,65	0,70	5,50	M2,0	20 Ncm



19.634	TIN screw hex	2,55	0°	7,85	1,65	0,70	5,50	M2,0	20 Ncm
--------	---------------	------	----	------	------	------	------	------	--------



INTERNAL HEX "CLICK"

compatible with
3i Certain®


SOCKET
HEX 1,20mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.011	Screw hex	2,00	0°	8,40	1,90	3,80	2,70	M1,6	20 Ncm



19.611	TIN screw hex	2,00	0°	8,40	1,90	3,80	2,70	M1,6	20 Ncm
--------	---------------	------	----	------	------	------	------	------	--------



A. HEAD DIAMETER (mm)

B. HEAD BASE ANGLE

C. OVERALL LENGTH (mm)

D. HEAD LENGTH (mm)

E. DISTANCE FROM
HEAD TO THREAD (mm)

F. THREAD LENGTH (mm)

G. THREAD KIND



INTERNAL HEX USA

compatible with
Zimmer Screw-Vent®

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.012	Screw hex	2,30	19°	8,20	2,55	2,98	2,60	UNF 1-72	30 Ncm
	19.042	Screw hex	1,95	19°	8,20	2,49	2,95	2,76	UNF 1-72	30 Ncm
	19.642	TIN screw hex	1,95	19°	8,20	2,49	2,95	2,76	UNF 1-72	30 Ncm



INTERNAL HEX MI

compatible with
MIS® Seven

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.056	Screw hex 3.0	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 Ncm
	19.081	Screw hex SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 Ncm
	19.656	TIN screw hex 3.0	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 Ncm
	19.681	TIN screw hex SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 Ncm
	19.110	Screw hex RP	2,00	0°	3,45	1,35	1,35	2,10	M1,4	15 Ncm



CONIC MC

compatible with
MIS® C1 internal

SOCKET
HEX 1,27mm




CONIC MV

compatible with
MIS® V3

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.056	Screw hex NP	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 Ncm
	19.081	Screw hex SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 Ncm
	19.656	TIN screw hex NP	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 Ncm
	19.681	TIN screw hex SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM
HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



EXTERNAL HEX BH

compatible with
BioHorizons® External

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.013	Screw hex NP	2,50	0°	7,00	1,00	2,10	3,90	UNF 1-72	30 Ncm
	19.014	Screw hex RP/WP	2,65	0°	7,10	1,00	2,10	4,00	M2,0	30 Ncm
	19.613	TIN screw hex NP	2,50	0°	7,00	1,00	2,10	3,90	UNF 1-72	30 Ncm
	19.614	TIN screw hex RP/WP	2,65	0°	7,10	1,00	2,10	4,00	M2,0	30 Ncm



INTERNAL HEX BH

compatible with
BioHorizons® internal

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.043	Screw hex	2,40	45°	7,24	2,42	2,36	2,45	UNF 1-72	30 Ncm
	19.643	TIN screw hex	2,40	45°	7,24	2,42	2,36	2,45	UNF 1-72	30 Ncm




INTERNAL CAM

compatible with
Camlog®

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.044	Screw hex NP	2,10	45°	10,45	2,35	4,60	3,50	M1,6	20 Ncm
	19.045	Screw hex WP	2,50	45°	10,45	2,35	4,60	3,50	M2,0	20 Ncm
	19.644	TIN screw hex NP	2,10	45°	10,45	2,35	4,60	3,50	M1,6	20 Ncm
	19.645	TIN screw hex WP	2,50	45°	10,45	2,35	4,60	3,50	M2,0	20 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



INTERNAL CON

compatible with
Conelog®

SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.074	Screw hex 3.3/3.8/4.3	2,26	30°	8,87	2,08	3,99	2,80	M1,6	20 Ncm
	19.076	Screw hex 5.0	2,67	30°	8,87	2,08	3,99	2,80	M2,0	20 Ncm
	19.674	TIN screw hex 3.3/3.8/4.3	2,26	30°	8,87	2,08	3,99	2,80	M1,6	20 Ncm
	19.676	TIN screw hex 5.0	2,67	30°	8,87	2,08	3,99	2,80	M2,0	20 Ncm



CONIC ANYR

compatible with
Megagen AnyRidge®

SOCKET
HEX 1,20mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.057	Screw hex	2,10	0°	8,00	2,00	2,78	3,22	M1,8	35 Ncm
	19.657	TIN screw hex	2,10	0°	8,00	2,00	2,78	3,22	M1,8	35 Ncm




CONIC AXIO

compatible with
Anthogyr AXIOM®

SOCKET
HEX 1,18mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.062	Screw hex Multi-Unit	2,18	0°	4,45	2,00	0,00	2,45	M1,4	15 Ncm
	19.077	Screw hex	2,22	70°	7,90	2,40	3,30	2,20	M1,6	25 Ncm
	19.662	TIN screw hex Multi-Unit	2,18	0°	4,45	2,00	0,00	2,45	M1,4	15 Ncm
	19.677	TIN screw hex	2,22	70°	7,90	2,40	3,30	2,20	M1,6	25 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



CONIC BIO

compatible with
Biotech KONTACT®

SOCKET
HEX 1,22mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.063	Screw hex NrW	2,03	30°	4,50	1,80	0,00	2,70	M1,6	20 Ncm
19.064	Screw hex Std	1,93	45°	4,50	1,95	0,00	2,55	M1,4	20 Ncm
19.078	Screw hex	1,93	45°	7,58	1,56	3,27	2,75	M1,4	20 Ncm
19.663	TIN screw hex NrW	2,03	30°	4,50	1,80	0,00	2,70	M1,6	20 Ncm
19.664	TIN screw hex Std	1,93	45°	4,50	1,95	0,00	2,55	M1,4	20 Ncm
19.678	TIN screw	1,93	45°	7,58	3,27	3,27	2,75	M1,4	20 Ncm



CONIC DENT

compatible with
Dentium Superline™ & Implantium®

SOCKET
HEX 1,28mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.066	Screw hex MUA	2,30	45°	5,00	1,93	0,57	2,50	M1,8	10 Ncm
19.090	Screw hex S-Line	2,30	33°	10,2	2,05	4,85	3,30	M2,0	30 Ncm
19.666	TIN screw hex MUA	2,30	45°	5,00	1,93	0,57	2,50	M1,8	10 Ncm
19.690	TIN screw hex S-Line	2,30	33°	10,2	2,05	4,85	3,30	M2,0	30 Ncm



CONIC IC

compatible with
Medentis ICX®

SOCKET
HEX 1,40mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.080	Screw hex	2,47	75°	9,10	2,97	4,13	2,00	M1,6	30 Ncm
19.680	TIN screw hex	2,47	75°	9,10	2,97	4,13	2,00	M1,6	30 Ncm




INTERNAL SEM

compatible with
Bego Semados®

SOCKET
HEX 1,27mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.091	Screw hex 3.25/5.5	2,41	55°	7,65	2,47	1,38	3,80	M1,8	30 Ncm
19.691	TIN screw hex 3.25/5.5	2,41	55°	7,65	2,47	1,38	3,80	M1,8	30 Ncm



A. HEAD DIAMETER (mm)
 B. HEAD BASE ANGLE
 C. OVERALL LENGTH (mm)
 D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
 F. THREAD LENGTH (mm)
 G. THREAD KIND



CONIC INK

compatible with
GlobalID IN-KONE®

SOCKET
HEX 1,22mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.035	Screw hex Straight	2,22	45°	5,00	2,10	0,00	2,90	M1,8	15 Ncm



19.036	Screw hex Angled	2,10	45°	4,80	2,10	0,00	2,70	M1,4	15 Ncm
--------	------------------	------	-----	------	------	------	------	------	--------



19.086	Screw hex	2,17	31°	8,45	2,24	4,21	2,00	M1,8	15 Ncm
--------	-----------	------	-----	------	------	------	------	------	--------



19.635	TIN screw hex Straight	2,22	45°	5,00	2,10	0,00	2,90	M1,8	15 Ncm
--------	------------------------	------	-----	------	------	------	------	------	--------



19.636	TIN screw hex Angled	2,10	45°	4,80	2,10	0,00	2,70	M1,4	15 Ncm
--------	----------------------	------	-----	------	------	------	------	------	--------



19.686	TIN screw	2,17	31°	8,45	2,24	4,21	2,00	M1,8	15 Ncm
--------	-----------	------	-----	------	------	------	------	------	--------



CONIC OSS

compatible with
Osstem® TS / HIOSSEN® ET

SOCKET
HEX 1,22mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.072	Screw hex Mini	2,22	45°	10,18	2,67	5,41	2,10	M1,6	20 Ncm



19.073	Screw hex Standard	2,33	45°	8,35	2,60	2,95	2,80	M2,0	30 Ncm
--------	--------------------	------	-----	------	------	------	------	------	--------



19.272	DLC screw hex Mini	2,22	45°	10,18	2,67	5,41	2,10	M1,6	20 Ncm
--------	--------------------	------	-----	-------	------	------	------	------	--------



19.273	DLC screw hex Standard	2,33	45°	8,35	2,60	2,95	2,80	M2,0	30 Ncm
--------	------------------------	------	-----	------	------	------	------	------	--------



19.672	TIN screw hex Mini	2,22	45°	10,18	2,67	5,41	2,10	M1,6	20 Ncm
--------	--------------------	------	-----	-------	------	------	------	------	--------



19.673	TIN screw hex Standard	2,33	45°	8,35	2,60	2,95	2,80	M2,0	30 Ncm
--------	------------------------	------	-----	------	------	------	------	------	--------



INTERNAL BT

compatible with
BTI®

SOCKET
HEX 1,22mm



REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
19.065	Screw hex MU	2,00	0°	3,62	1,50	0,32	1,80	M1,4	10 Ncm




19.079	Screw hex	2,55	0°	7,05	1,80	2,05	3,20	M1,8	35 Ncm
--------	-----------	------	----	------	------	------	------	------	--------



19.665	TIN screw	2,00	0°	3,62	1,50	0,32	1,80	M1,4	10 Ncm
--------	-----------	------	----	------	------	------	------	------	--------



19.679	TIN screw MU	2,55	30°	7,05	1,80	2,05	3,20	M1,8	35 Ncm
--------	--------------	------	-----	------	------	------	------	------	--------



A. HEAD DIAMETER (mm)

B. HEAD BASE ANGLE

C. OVERALL LENGTH (mm)

D. HEAD LENGTH (mm)

E. DISTANCE FROM
HEAD TO THREAD (mm)

F. THREAD LENGTH (mm)

G. THREAD KIND



MULTI-UNIT CONNECTION

SOCKET
HEX 1,20mm



SOCKET
HEX 1,27mm



SOCKET
UniGrip®



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.006	Screw UniGrip® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.031	Screw UniGrip® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 Ncm
	19.055	Screw hex 1,20	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.065	Screw hex 1,22	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.085	Screw hex Multi-Unit 4.8	2,15	0°	3,50	1,30	0,00	2,20	M1,6	25 Ncm
	19.110	Screw hex 1,27 RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.206	DLC screw UniGrip® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.231	DLC screw UniGrip® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 Ncm
	19.606	TIN screw hex 1,20 on 3i® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.631	TIN screw UniGrip® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 Ncm
	19.685	TIN screw hex Multi-Unit 4.8	2,15	0°	3,50	1,30	0,00	2,20	M1,6	25 Ncm

Screw C-Base®


SOCKET
HEX 1,00mm



SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.052	Screw hex 1,00 C-Base®	2,15	45°	8,87	1,91	4,16	2,80	M1,6	20 Ncm
	19.702	Screw hex 1,27 C-Base®	1,90	60°	9,37	1,50	5,97	3,40	M1,8	15 Ncm



A. HEAD DIAMETER (mm)
B. HEAD BASE ANGLE
C. OVERALL LENGTH (mm)
D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)
F. THREAD LENGTH (mm)
G. THREAD KIND



MULTI-UNIT CONNECTION

SOCKET
HEX 1,20mm



SOCKET
HEX 1,22mm



SOCKET
HEX 1,27mm



SOCKET
UniGrip®



SOCKET
Torx



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.006	Screw UniGrip® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.031	Screw UniGrip® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 Ncm
	19.055	Screw hex 1,20	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.065	Screw hex 1,22	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.085	Screw hex Multi-Unit 4.8	2,15	0°	3,50	1,30	0,00	2,20	M1,6	15 Ncm
	19.110	Screw hex 1,27 RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.206	DLC screw UniGrip® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.231	DLC screw UniGrip® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 Ncm
	19.606	TIN screw hex 1,20 on 3i® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 Ncm
	19.631	TIN screw UniGrip® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 Ncm
	19.685	TIN screw hex Multi-Unit 4.8	2,15	0°	3,50	1,30	0,00	2,20	M1,6	15 Ncm
	19.069	AOX Screw Torx	2,20	0°	4,10	1,40	0,35	2,35	M1,4	15 Ncm
	19.093	AOX Screw hex 1,20	2,20	0°	4,10	1,40	0,35	2,35	M1,4	15 Ncm
	19.094	AOX Screw hex 1,22	2,20	0°	4,10	1,40	0,35	2,35	M1,4	15 Ncm
	19.095	AOX Screw hex 1,27	2,20	0°	4,10	1,40	0,35	2,35	M1,4	15 Ncm
	19.096	AOX Screw UG	2,20	0°	4,10	1,40	0,35	2,35	M1,4	15 Ncm
	19.022	AOX Screw NEO	2,20	0°	4,10	1,40	0,35	2,35	M1,4	15 Ncm
	19.269	AOX Screw Torx DLC	2,20	0°	4,10	1,40	0,35	2,35	M1,4	15 Ncm

Screw C-Base®


SOCKET
HEX 1,00mm



SOCKET
HEX 1,27mm



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.052	Screw hex 1,00 C-Base®	2,15	45°	8,87	1,91	4,16	2,80	M1,6	20 Ncm
	19.702	Screw hex 1,27 C-Base®	1,90	60°	9,37	1,50	5,97	3,40	M1,8	15 Ncm



- A. HEAD DIAMETER (mm)
- B. HEAD BASE ANGLE
- C. OVERALL LENGTH (mm)
- D. HEAD LENGTH (mm)


- E. DISTANCE FROM HEAD TO THREAD (mm)
- F. THREAD LENGTH (mm)
- G. THREAD KIND

SCREW ANGLEBase®

SOCKET
Torx



	REF	DESCRIPTION	A	B	C	D	E	F	G	TORQUE
	19.306	Screw MUA ANGLEBase®	2,20	0°	3,60	1,50	0,10	2,00	M1,4	15 Ncm
	19.331	Screw UNI ANGLEBase®	2,30	0°	4,10	1,40	0,70	2,00	M1,8	15 Ncm
	19.400	Screw ANGLEBase®	2,30	0°	7,40	1,40	3,30	2,70	UNF1,72	30 Ncm
	19.401	Screw ANGLEBase®	2,10	0°	6,10	1,40	2,50	2,20	M1,6	35 Ncm
	19.402	Screw ANGLEBase®	2,50	0°	5,90	1,40	1,69	2,81	M2,0	35 Ncm
	19.403	Screw ANGLEBase®	2,30	0°	5,90	1,40	2,50	2,00	M1,8	35 Ncm
	19.404	Screw ANGLEBase®	2,30	0°	7,20	1,40	3,55	2,25	M1,8	35 Ncm
	19.405	Screw ANGLEBase®	2,45	0°	8,05	1,40	4,60	2,05	M2,0	35 Ncm
	19.417	Screw ANGLEBase®	2,30	0°	4,90	1,40	2,90	3,85	UNF-1-72	30 Ncm
	19.425	Screw ANGLEBase®	2,45	0°	8,90	1,40	5,45	2,00	M2,0	25 Ncm
	19.438	Screw ANGLEBase®	2,10	0°	7,90	1,40	4,50	2,05	M1,6	● Ncm
	19.441	Screw ANGLEBase®	2,10	0°	7,30	1,40	3,70	2,00	M1,6	35 Ncm
	19.442	Screw ANGLEBase®	2,45	0°	7,20	1,40	3,75	2,05	M2,0	35 Ncm
	19.443	Screw ANGLEBase®	2,10	0°	8,15	1,40	4,75	2,00	M1,6	35 Ncm
	19.444	Screw ANGLEBase®	2,10	0°	9,20	1,40	5,60	2,20	M1,6	● Ncm
	19.446	Screw ANGLEBase®	2,45	0°	6,70	1,40	3,25	2,05	M2,0	35 Ncm
	19.460	Screw ANGLEBase®	2,30	0°	7,42	1,40	4,02	2,00	M1,8	25 Ncm
	19.461	Screw ANGLEBase®	2,45	0°	7,75	1,40	4,30	2,05	M2,0	25 Ncm
	19.468	Screw ANGLEBase®	2,10	0°	6,90	1,40	3,05	2,45	M1,4	15 Ncm
	19.474	Screw ANGLEBase®	2,10	0°	8,70	1,40	5,30	2,00	M1,6	20 Ncm
	19.476	Screw ANGLEBase®	1,90	0°	8,50	1,40	5,05	2,05	M2,0	20 Ncm



● Internal Hex "Click"= 25 Ncm
Internal Hex FD= 24 Ncm

● Internal CAM= 20 Ncm
Internal OSS (Mini)= 24 Ncm
Internal ICX= 30 Ncm

A. HEAD DIAMETER (mm)

B. HEAD BASE ANGLE

C. OVERALL LENGTH (mm)

D. HEAD LENGTH (mm)

E. DISTANCE FROM HEAD TO THREAD (mm)

F. THREAD LENGTH (mm)

G. THREAD KIND

